LLC RC ART
Ekaterinburg, Russia

Transcutaneous Electrostimulator

Operations Manual

Osteo · DENS

RC ART 08.0-03.71-01 RE

LLC „RC ART“
Ekaterinburg, Russland

Elektrostimulator

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Part 2. User's Instructions

1. Main Performance Data.....	30
2. Function.....	32
3. Indications and Contraindications to Application.....	33
4. Conditions of Treatment.....	34
5. Intensity of the Electrostimulating Treatment.....	35
6. Recommendations for Application.....	37

Coupon for Warranty Repair.....	81
--	-----------


Certificate of Acceptance.....	83
---------------------------------------	-----------

This Operations Manual is intended for the Osteo-DENS electrostimulator with built-in electrodes for stimulation of reflex zones on the posterior surface of the human neck.

The Operations Manual includes the Technical Passport (part 1) and User's Instructions (part 2).

PART 1

Technical passport

1. SAFETY

All information in this operations manual contains important information about your safety and recommendations on proper use and maintenance of the device.

The device does not constitute any danger because of internal low voltage power source isolated from the body of the apparatus (article of type B and body of type F).

The apparatus must not be used for treatment of patients with implanted electronic devices (for example, pacemaker) and for treatment of patients who have individual electric current intolerance.

Use of the device in direct front projection of the heart is prohibited.

Don't switch the patient to any high-frequency electric device during stimulation. Simultaneous use of the device and other electric equipment can result in burns and possible damage of the device.

Work near short-wave and microwave equipment can result in instability of output parameters of the device.

Avoid long-term exposure of the device to direct sunlight at high (>25°) air temperature. This can result in malfunction of the device indicator or damage of the power elements.


The device has fragile elements. Keep it safe from blows.


The device is not waterproof. Keep it safe from water penetration.


All repair works on the device must be provided by qualified specialists at the manufacturer.

Transportation conditions: temperature -50°C to +50°C, relative air humidity 30 to 93%, atmospheric pressure 70 to 106 kPa.

Storage conditions: temperature -50°C to +40°C, relative air humidity 30 to 93%, atmospheric pressure 70 to 106 kPa.

Operation conditions: temperature +10°C to +35°C, relative air humidity 30 to 93%, atmospheric pressure 70 to 106 kPa.

Attention! *In case of storing the device at temperature below 10°C, keep it under normal conditions minimum two hours before use.*

Recycling: All packaging materials are environmental-friendly and can be recycled.

Individual assembly of electrical and electronic equipment.

The device contains valuable materials that may be used again after recycling considering the requirements for environment preservation. Deliver it to special service centers (consult with appropriate services of your region) for collection and recycling.

2. INTENDED USE

The Osteo-DENS device is intended for electrostimulation of reflex zones on the posterior surface of the human neck. The device has built-in electrodes.

The Osteo-DENS apparatus is intended for individual application in patient-care institutions and in life conditions in compliance with the recommendations of the attending doctor.

3. SPECIFICATIONS

3.1. Electric impulse parameters under nominal power voltage have the following output parameters:

3.1.1. Under minimum signal power with load of 2 kOhm the duration of the positive impulse is $25 \pm 10 \mu\text{s}$.

Amplitude of the positive impulse is maximum 30 V.

Amplitude of the negative impulse is maximum 30 V.

3.1.2. Under maximal signal power with load of 20 kOhm the duration of the positive impulse $250 \pm 70 \mu\text{s}$
Amplitude of the positive impulse is $30 \pm 10 \text{ V}$.










Amplitude of the negative impulse with load ($20 \pm 1 \text{ kOhm}$) – $150 \pm 70 \text{ V}$.


3.1.3. Amplitude of treatment at the minimum power level is $\approx 50\%$ of the amplitude at the maximum power level ($R = 20 \text{ kOhm}$).

3.2. Dependence of the impulse shape on the load resistance at the medium and maximum power

Table 1

Load resistance	Power level 7 units		Power level 15 units	
	Shape of the output signal	V_{P-P}	Shape of the output signal	V_{P-P}
No load		$\approx 140 \text{ V} \pm 20\%$		$\approx 140 \text{ V} \pm 20\%$

<p>200 Ohm</p>		<p>$\approx 9\text{ V}$ $\pm 20\%$</p>		<p>$\approx 10\text{ V}$ $\pm 20\%$</p>
<p>500 Ohm</p>		<p>$\approx 24\text{ V}$ $\pm 20\%$ leff \approx 16 ns $E_{\text{pulse}} \approx$ 42 mJ</p>		<p>$\approx 28\text{ V}$ $\pm 20\%$ leff \approx 8 ns $E_{\text{pulse}} \approx$ 78 mJ</p>
<p>1 kOhm</p>		<p>$\approx 42\text{ V}$ $\pm 20\%$</p>		<p>$\approx 54\text{ V}$ $\pm 20\%$</p>
<p>2 kOhm</p>		<p>$\approx 60\text{ V}$ $\pm 20\%$</p>		<p>$\approx 84\text{ V}$ $\pm 20\%$</p>
<p>10 kOhm</p>		<p>$\approx 140\text{ V}$ $\pm 20\%$</p>		<p>$\approx 140\text{ V}$ $\pm 20\%$</p>

Notes: curves in phase A of mode 1 at the point of maximum power.

3.3. The apparatus provides for selection and setting of two electrostimulation modes.

3.4. Maximum consumable current (with 3 V voltage) in the device is maximum 140 mA.

3.5. Power source: 2 batteries type LR6/AA, total voltage $3.0 \pm 0.6\text{ V}$.

It is allowed to use accumulators of the similar type, size and voltage*.

3.6. Weight of the main parts of the electrostimulator:
— control unit, maximum 80 gr.
— block of neck electrodes, maximum 170 gr.

3.7. Dimensions of the main parts of the electrostimulator:
— control unit, maximum 130x50x22 mm;
— block of neck electrodes, maximum 120x60x40 mm.

* *order of operation (types of chargers, charging methods) is provided in the manual for accumulators. Period of apparatus operation with accumulators depends on the accumulators' specifications.*

3.8. The device will automatically switch off not later than in 2 minutes after the last pressing of the apparatus buttons (except for On/Off button) or last application of the electrodes to the patient's skin surface.

3.9. Electromagnetic Radiation

Table 2

Test	Compliance	Application Conditions
HF radiation CISPR 11	Class B	Osteo-DENS is suitable for use in all establishments including home application.

3.10. HF Resistance

Table 3

Test	IEC 60601-1-2 test conditions	Permissible level
IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms
IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m

3.11. Resistance to Electromagnetic Fields

Table 4

Test	Test Level	Level of Compliance	Application Conditions
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV infl.	±4kV contact ±8kV infl.	Floors to be wooden, concrete or ceramic tile. If floors are covered with synthetic material, relative air humidity to be minimum 40%.

Magnetic fields IEC 61000-4-8	3A/m	3A/m	Magnetic fields parameters to be typical, for business buildings and medical and prophylactic institutions.
----------------------------------	------	------	---

3.12. Recommendations on determining the required distance between the Osteo-DENS electrostimulator and radio equipment

Table 5

Stated maximum output power of the transmitter R (W)	Radiation Frequency of the Transmitter and Equation to Determine the Distance d (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	150 kHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

4. COMPLETE SET

The complete set of delivery of the Osteo-DENS device should correspond to Table 6.

Table 6

Name of the device, part, document	Quantity
Control unit	1
Block of neck electrodes	1
Operations Manual	1
Case	1
Power source (AA type)	2

5. DESIGN OF THE DEVICE

5.1. General Information

The Osteo-DENS apparatus consists of the block of neck electrodes and the control unit connected with the cable. The power sources to be placed into the compartment of the block of neck electrodes.

5.2. Block of neck electrodes

The block of neck electrodes (Fig. 1) is to be fixed on


Fig. 1. Appearance of the Block of Neck Electrodes

1. Body
2. Built-in Electrodes
3. Cover of the battery compartment

the posterior surface of the patient's neck with the help of a soft elastic bracing. Electrostimulation of the reflex zones on the posterior surface of the human neck is conducted directly from the block of neck electrodes.

Attention! *Before the treatment procedure with the apparatus it is recommended to take off current-conducting decorations (metallic chains and so on).*

Attention! *It is strictly forbidden to connect to the block of neck electrodes any devices except for the control unit in the complete set of the Osteo-DENS apparatus. Connection of an inappropriate device can result in malfunction of this device or in damage of neck electrodes and be dangerous for the user. Breakdown of neck electrodes resulting from connection of inappropriate devices shall not be a warranty event!*

Attention! *Allergic reactions are possible due to individual intolerance of the materials used in the apparatus. In case of allergic reactions you should stop the procedure and address an expert.*

5.3. Control Unit

The control unit (Fig. 2) provides control of modes and power of electrostimulation in the block of neck electrodes via the connection cable.


Fig. 2. Appearance of the control unit

Attention! *The control unit can be connected to the block of neck electrodes in the complete set of the Osteo-DENS apparatus. Connection of the control unit to an inappropriate device can result in malfunction of this device or in damage of the control unit. Breakdown of the control unit resulting from connection of inappropriate devices shall not be a warranty event!*


5.3.1. Control Buttons on the keyboard

Table 7

Control Buttons*			
			
Switch on/off of the apparatus	Switching between electrostimulation modes	Reduction of the treatment power	Increase of the treatment power

5.3.2. Indicator of the Control Unit

1. Digital indicator and scale of duration of the therapeutic signal influence

* When pressing any of the keyboard buttons, the indicator flashes.

6. OPERATION MODE


6.1. Preparing the apparatus for operation

Open the cover of the battery compartment on the block of neck electrodes (BNE). Place the batteries as per item 8 of this Technical Passport. Close the battery compartment.

Connect the BNE and control unit (CU) with the cable. Fix the BNE on your posterior neck surface yourself or ask for help, electrodes should fit on the skin but the bracing shall not be too tight. For better contact the skin in the area to be treated can be wiped with a tampon wetted with water or apply some Malavtilin cream (the procedure may be started after full absorption of the cream).

6.2. Switching the apparatus on

Shortly press the  button on the control unit, you will hear the sound signal with introduction and the indicator will show the following information:

6.3. Application of the apparatus

6.3.1. To change the therapy mode shortly press the  button.

To set or change the electrostimulation power use  and  buttons.

Attention! *The power increase is controlled subjectively following the patient's sensations upon contact of the electrodes with the skin surface. Do not surpass the pain threshold.*

At the power level more than zero and with the contact of all electrodes with the skin surface, the apparatus will switch to the mode of electrostimulation in compliance with the selected mode and power level.

If in the process of therapy at least one of the electrodes loses contact with the skin surface, you will hear a sound signal corresponding to this event and the procedure will stop for 3 seconds. If during this time contact with skin is restored, the procedure will resume; in case the contact is not restored, the apparatus will switch into the stand-by mode.

During the procedure it is allowed changing the intensity of the electric stimulating impulse and the procedure is neither paused nor resumed.

The procedure will automatically stop upon expiration of the treatment time for the mode applied. When the procedure stops, the actual power of electrostimulation reduces to the minimum value, but the indicator will show the power level set before.

You can repeat the procedure either by changing the mode (one or two times) with the  button, or by short taking off and again applying electrode to the skin, or by changing the treatment power.

For forced switch off of the therapy mode, either reduce the treatment power to zero or remove the block of neck electrodes from the patient's neck.

6.4. Switching the apparatus off

To switch the apparatus off press and hold (more than 3 seconds) the  button until you hear the signal indicating the end of the procedure, after that release the  button and after the end of the sound signal the apparatus will switch off.

The apparatus shall switch off automatically if it is not applied during 2 minutes. The switched on apparatus is considered not to be use if there is no electrodes contact with the skin surface and none of the buttons is pressed.

7. TECHNICAL MAINTENANCE

7.1. Daily technical maintenance should include the following:

- external examination of the apparatus;
- disinfection of electrodes.

Use standard disinfection means (e.g. 70% alcoholic solution) and soft napless napkins to clean the electrodes.


7.2. If the apparatus is supposed not to be used for a long period (more than 3 days), remove the power sources from the compartment of the block of neck electrodes (BNE).

7.3. Having seen the blinking battery symbol, change the power sources.

8. CHANGE OF POWER SOURCES

Change of power sources:

- open the battery compartment;
- take the power source out;
- let the apparatus stay without power sources minimum 2 minutes;
- set new power sources*, follow the polarity.

-
- * Set only those power sources that are provided for this device – type LR6/AA, voltage rating 1.5 V, or appropriate accumulators with nominal voltage 1.2 V.

9. TROUBLESHOOTING

Possible problems and troubleshooting are provided in Table 8.

Table 8

Trouble	Possible reason	Method of eradication
The apparatus does not switch on after pressing the On/Off button	Power sources are missing	Set a new set of power sources (see section "Change of Power Sources")
	The polarity is mixed when power sources setting into the compartment	Check and, if necessary, set power sources per marking on the apparatus body
	Voltage of power sources is less than 2.2 V	Replace power sources (see section "Change of Power Sources")
The device produces sound signals during switching on and automatically switches off	Voltage of power sources is less than 2.2 V	Replace power sources (see section "Change of Power Sources")
The apparatus does not switch to the Therapy mode	The contact of electrodes with the skin surface is not enough	Fix the apparatus electrodes on the skin surface. If necessary, wet the skin surface in the area of contact with electrodes

<p>The apparatus does not switch off automatically when there is no contact of electrodes with the skin surface and not touching the controls more than 2 minutes</p>	<p>Current leakage on the therapeutic electrodes of the apparatus</p>	<p>Clean and wipe the apparatus therapeutic electrodes with an alcohol or alcohol-containing liquid</p>
<p>The apparatus switches off itself during the procedure</p>	<p>Voltage of power sources is less than 2.2 V</p>	<p>Replace power sources (see section "Change of Power Sources")</p>
	<p>Bad contact of the cable in the plug of the block of neck electrodes</p>	<p>Check that the control unit is reliably connected to the block of neck electrodes</p>
<p>Power sources discharge very quickly</p>	<p>Bad quality of power sources</p>	<p>Use power sources of good quality (we recommend alkaline power sources) or accumulator of proper type and size of voltage maximum 1.5 V</p>

10. WARRANTY CONDITIONS

10.1. The operation lifetime is 5 years.

Observation of operation conditions can considerably increase the lifetime set by the manufacturer officially.

10.2. Warranty period of operation is 24 months from the date of sale.

10.3. The seller (manufacturer) or organization carrying out functions of the seller (manufacturer) on a contractual basis is not responsible for defects should they occur after the disposal of the device as a result of:

- 1) a failure on the part of the consumer to comply with rules of transportation, storage, maintenance and operation provided in the present manual;
- 2) mechanical damages;
- 3) actions of the third party;
- 4) force-majeure.

10.4. Warranty obligations do not apply to products with broken manufacturer's seals.

10.5. In case of unit breakdown or malfunction within the warranty period, as well as in case of incomplete shipping is found, the owner must send the following documents to the manufacturer's address or manufacturers' representative: claim for repair (exchange) with

name, address, telephone number; defects list with brief description of the malfunction, date and conditions of its appearance.

Manufacturer's address:

LLC "RC ART", 620146 Russia
Ekaterinburg, Akademika Postovskogo Str., 15
Phone: +7 (343) 267-23-30
<http://www.denascorp.ru>,
e-mail: corp@denascorp.ru

Official Representative in the European Union:

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
Trade representative in the European Union:

DENAS-CZ s.r.o., Czech Republik, 360 01, Karlovy Vary,
SHOPPING CENTER "atrium", Karla IV. 505/1, office 209,
phone: (+420) 353 549 285; fax: (+420) 359 019 209

PART 2

User's instruction

1. MAIN PERFORMANCE DATA

The neck is the most mobile and vulnerable part of the spine that is why every person has pain in the neck during his or her life regardless of the age and gender. Among other reasons for pain the most frequent ones are diseases and traumas of the cervical spine, muscles, and neck ligaments.

The name of the cranial vertebra is atlas. Atlas in the Greek mythology is a Titan who supports the dome of heaven on his shoulders. The neck supports the head, and blood supply of the brain and general state of health depend on health of the cervical spine. That is why, pain in the neck needs special attention and special treatment.

Among methods of treatment of acute and chronic pain in the neck, non-drug methods are preferred nowadays. No matter which reason results in pain in the neck, a complex of proper actions is needed for efficient treatment and prophylaxis of attacks, among which dynamic electroneurostimulation (DENS) is holding an important position.

Many researches evidence that DENS treatment effect is based on multilevel reflex and neurochemical reactions which activate a chain of regulatory and adaptive

mechanisms of the body. As a result, pain syndromes are eliminated, blood circulation is improved, biologically active substances formation and metabolic processes in tissues are activated, muscle and vascular tone are normalized. The dynamic electroneurostimulation promotes improvement of the general state of health, mood and increase of work efficiency.

2. FUNCTION

The apparatus is intended for correction of functional disorders conditioned by affection of the cervical spine, neck and shoulder girdle muscles.

3. INDICATIONS AND CONTRAINDICATIONS FOR APPLICATION

3.1. Indications for Application

Pain syndrome, muscle tension, restriction of movement, muscle fatigue and other functional disorders, conditioned by diseases of the cervical spine, neck muscles and shoulder girdle, and also long forced position or heavy exercises.

3.2. Contraindication for Application:

- individual intolerance of electric current;
- presence of an implanted pacemaker;
- epilepsy;
- neoplasms of any etiology and localization;
- acute fevers of unknown etiology;
- venous thrombosis;
- condition of acute mental excitement, alcohol or drug intoxication.

Attention! *In the above cases of contraindications application of the apparatus should be first discussed with the attending physician.*

Attention! *In case of individual intolerance when during the treatment procedure the state of health is worsened and pain increases, stop the procedure and urgently consult a doctor.*

4. CONDITIONS OF TREATMENT

You do not need any special conditions for the procedure. You can conduct the treatment procedures yourself. The therapy is taken in a comfortable sitting or lying position. After the treatment procedure the patient should relax for 10–15 minutes.

Attention! *After each procedure all electrodes to be treated with a standard disinfecting solution (e.g. 70% alcoholic solution). The device should be kept with dry electrodes.*

5. INTENSITY OF THE ELECTROSTIMULATING TREATMENT

The intensity (power) of electrostimulation is determined by subjective patient's sensations. Electrostimulation intensity is conditionally divided into three energy levels (electric ranges).

The first minimum level – the patient does not have any subjective sensations or feels subtle vibrations in the subelectrode zone. It is applied in cases when treatment should not be intensive, i.e. for patients with diseases of the cardio-vascular system (arterial hypertension, arterial hypotension, syndrome of vegetative-vascular dystonia). For patients who often have headaches, dizziness, treatment also shall be conducted at the minimum power level.

The second comfortable level – the patient feels light pricking, vibration or light burning without pain. It is applied for mild pains and medium intensity pains. This is the most frequently used power level.

The third maximum level – the patient feels painful pricking or burning. Such intensity can be accompanied by involuntary contraction of muscles close to the electrodes. It is used in case of expressed pain syndrome.

Attention! *Power of electrostimulation can be increased or reduced during treatment stages depending on the extent of the patient's sensitivity and as the pain syndrome reduces.*


Attention! *The power level is controlled subjectively following the patient's sensations. Do not surpass the pain threshold.*

6. RECOMMENDATIONS FOR APPLICATION


Attention! For patients with diseases of the cardio-vascular system (arterial hypertension, arterial hypotension, syndrome of vegetative-vascular dystonia) and for patients who often have headaches, it is recommended to conduct the procedure at the minimum power level only.

6.1. — A program (pain therapy).


The program is applied in cases when the main complaint is severe pain in the neck or shoulder girdle area.

For the procedure, it is necessary to:
1. Switch the apparatus on with the  button on the control unit.


2. Fix the block of neck electrodes in the area of posterior neck surface in such a way that both electrodes fit on the skin.

And on the control unit display you will see the indicators of the electrodes contact with skin. With good contact both indicators will switch on. If some of the electrodes does not fit on the skin, the related indicator will not switch on.

Attention! *In case of no contact of both electrodes with the skin during two minutes, the apparatus will switch off.*

3. By pressing the  button, select "Pain Therapy" program, you will see the  symbol on the display.

4. Set the required treatment power with  and  buttons of the control unit.

Attention! *The power level should be set on more than zero to start the program.*

5. Conduct the procedure.
The duration of this treatment program is 20 minutes. Upon the end of the treatment program the procedure will stop automatically.


Attention! *If you want to make a pause during the procedure, reduce the treatment power to zero.*

Repeat the procedure as needed in case when there are complaints.

Attention! *Keep the apparatus with dry electrodes.*


6.2.  — B program (relaxation).

The program is applied in cases when there is no pain syndrome or it is very mild, but the patient has restraint of neck movement, woodiness, fatigue, excessive strain of the shoulder girdle muscles or posterior neck surface.


For the procedure, it is necessary to:

1. Switch the apparatus on with the  button on the control unit.


2. Fix the block of neck electrodes in the area of posterior neck surface in such a way that both electrodes fit on the skin.

And on the control unit display you will see the indicators of the electrodes contact with skin. With good contact both indicators will switch on. If some of the electrodes does not fit on the skin, the related indicator will not switch on.

Attention! *In case of no contact of both electrodes with the skin during two minutes, the apparatus will switch off.*

3. By pressing the  button, select "Relaxation" program, you will see the  symbol on the display.

4. Set the required treatment power with  and  buttons of the control unit.

Attention! *The power level should be set on more than zero to start the program.*

5. Conduct the procedure.

The duration of this treatment program is 15 minutes. Upon the end of the treatment program the procedure will stop automatically.

Repeat the procedure as needed in case when there are complaints.

Attention! *Keep the apparatus with dry electrodes.*

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Date of sending for repair

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Reason for repair / Ursache Abgabe zur Reparatur _____

Note about repair / Reparaturvermerk _____

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Date of acceptance
Annahme am _____

Warranty for a repaired product is valid for 6 months beginning as of the moment when product is accepted following repairs. If a warranty period beginning on the date of purchasing is more than 6 months, warranty is recalculated for a longer period. Warranty period is increased accordingly by the time a product has spent at a repair station.

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CERTIFICATE OF ACCEPTANCE ANNAMEBESCHEINIGUNG

Electrostimulator with built-in electrodes for stimulation of biologically active points and zones Osteo-DENS

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I received information about warranty conditions, the apparatus was checked, and I have no claims to the complete set, external appearance of the apparatus.

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Signature of the customer

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Carefully examine the apparatus when buying!

Body and display defects (scratches, cracks, fractures) are not warranty events. Apparatuses with such defects are not exchangeable, beyond repair and not returnable.

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