

















Comparative table of DENAS/DiaDENS devices

Properties/functions/specifications	Devices for Personal Usage and for Doctors						Devices for Doctors only	
	DiaDENS-PCM	DiaDENS-T	DENAS	DiaDENS-Cardio	DENAS+	Osteo DENS	DiaDENS-DT	DiaDENS-PC
Appearance								
Purpose of use	Universal device for different health related problems. Fast pain relief, treatment of acute and especially chronic diseases	Universal device for different health related problems. Fast pain relief, treatment of acute and chronic diseases	Universal device for different health related problems. Pain relief and treatment of acute diseases	Specific device for treatment of hypertension (high blood pressure)	Universal device for different health related problems. Pain relief and treatment of acute diseases.	Specific device for treatment neck pain and shoulder girdle.	Universal device for different health related problems. It includes acupuncture diagnostic. Fast pain relief, treatment of acute and chronic diseases	Universal device for different health related problems. It includes acupuncture diagnostic. Fast pain relief, treatment of acute and especially chronic diseases
Ability to hook up special accessories for therapy (DENS-Glasses, DENS Applicator, Facial Probe, Hair Comb Probe, Therapeutic Acupuncture Electrode)	+ You need small adapter to hook up the following: DENS-Glasses, DENS Applicator, Therapeutic Acupuncture Electrode.	+ You need small adapter to hook up the following: Facial Probe and Hair Comb Probe.	+ You need adapter (it looks like cradle) to hook up all accessories.	-	+ You need small adapter to hook up the following: DENS-Glasses, DENS Applicator, Therapeutic Acupuncture Electrode.	-	+ You need small adapter to hook up the following: Facial Probe and Hair Comb Probe.	+ You need small adapter to hook up the following: DENS-Glasses, DENS Applicator, Therapeutic Acupuncture Electrode.

Ability to hook up special accessories for acupuncture diagnostic (acupuncture diagnostic probes: active and passive)	-	-	-	-	-	-	+	+
Ability to connect to PC for acupuncture diagnostic	-	-	-	-	-	-	-	+
LCD Screen	+ LCD Screen with backlight	+	-	+	-	+ LCD Screen with backlight	+	+
Energy Level Control	Digital Control with levels from 0 up to 99	Digital Control with levels from 0 up to 99	LED ruler with 10 Energy Levels.	- Fixed Energy Level	LED ruler with 10 Energy Levels.	Digital Control with levels from 0 up to 15	Digital Control with levels from 0 up to 99	Digital Control with levels from 0 up to 99
Form of build-in electrode (on-board electrode)								
Frequencies	Standard frequencies: 10, 20, 60, 77, 140, 200 Hz Frequency Modulation: «77 10» Hz (special regimen with alternate frequencies 77	Standard frequencies: 10, 20, 60, 77, 140, 200 Hz	Standard frequencies: 10, 77 Hz	One ultra low frequency: 9.2 Hz and Frquency Modulation: «77 10» Hz	Standard frequencies: 10, 77 Hz	-	Standard frequencies: 10, 20, 60, 77, 140, 200 Hz	Standard frequencies: 10, 20, 60, 77, 140, 200 Hz Frequency Modulation: «77 10» Hz (special regimen with alternate frequencies 77

	<p>Hz and 10 Hz) Amplitude Modulation: «77AM» (special regimen with variable amplitude) Ultra low frequencies: 1-9.9 Hz with step 0.1 Hz.</p>							<p>Hz and 10 Hz) Amplitude Modulation: «77AM» (special regimen with variable amplitude) Ultra low frequencies: 1-9.9 Hz with step 0.1 Hz.</p>
<p>Regimens of work:</p>	<p>1. «Therapy» + part of this regimen - «Express» for treating 15 health related problems 2. «MED» for prophylaxis and treating General Zones 3. «Screening» regimen for diagnostic pathological spots.</p>	<p>1. «Therapy» 2. «MED» for prophylaxis and treating General Zones 3. «Screening» regimen for diagnostic pathological spots.</p>					<p>1. «Therapy» 2. «MED» for prophylaxis and treating General Zones 3. «Screening» regimen for diagnostic pathological spots. 4. «Voll» - for acupuncture diagnostic 5. «Bioreper» - acupuncture diagnostic for ears</p>	<p>1. «Therapy» + part of this regimen - «Express» for treating 15 health related problems 2. «MED» for prophylaxis and treating General Zones 3. «Screening» regimen for diagnostic pathological spots. 4. «Voll» - for acupuncture diagnostic 5. «BioVoll» - advanced Voll diagnostic 6. «Bioreper» - acupuncture diagnostic for ears 7. «MiniAS» —</p>

								acupuncture diagnostic for Su Jok projections
Complete set	<ol style="list-style-type: none"> 1. DiaDENS-PCM device 2. Operations manual including a technical passport and user’s instructions 3. Package type AAx2 5. Holder lace 	<ol style="list-style-type: none"> 1. DiaDENS-T device 2. Operation manual 3. Therapeutic remote electrode, (acupuncture probe) 4. Container 5. Package 6. Power source type 6F22 	<ol style="list-style-type: none"> 1. DENAS device 2. Operations manual including a technical passport and user’s instructions (very brief) 3. Container 4. Package 5. Power source type AA 	<ol style="list-style-type: none"> 1. DiaDENS Cardio device 2. Operations manual including a technical passport and user’s instructions 3. Plastic Container 4. Package 5. Power source type AAx2 	<ol style="list-style-type: none"> 1. DENAS+ device 2. Operations manual including a technical passport and user’s instructions (very brief) 3. Soft case 4. Package 5. Power source type AAx2 6. Adapter (for connection to some accessories) 	<ol style="list-style-type: none"> 1. OsteoDENS device 2. Remote Controller with LCD 3. Operations manual including a technical passport and user’s instructions 4. Four Small pillows 5. Package 6. Power source type AAx2 	<ol style="list-style-type: none"> 1. DiaDENS-DT device 2. Operation manual 3. Therapeutic remote electrode, (acupuncture probe) 4. Diagnostic remote electrode (acupuncture probe) 5. Container 6. Package 7. Power source type 6F22 	<ol style="list-style-type: none"> 1. “DiaDENS-PC” device 2. Operation Manual, including Technical Passport and Operation Instructions 3. Therapeutic remote electrode, (acupuncture probe) 4. Diagnostic remote electrode (acupuncture diagnostic probe) 5. Connecting computer cable 6. CD, software 7. User Manual for the DiaDENS-PC software 8. Case 9. Cover for device 10. Packaging 11. Power source type 6F22